

# HOXD12 antibody



Catalog Number: 111438

## Product name

HOXD12 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

HOXD12 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, L02 cells, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 30kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: BC121103). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111438(HOXD12 antibody) at dilution of 1:200