

HOXD13 antibody



Catalog Number: 111439

Product name

HOXD13 antibody

Specificity

Human; other species not tested.

Antibody description

HOXD13 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, HL-60 cells. Positive IF detected in A549 cells. Observed molecular weight by Western-blot: 38-40 kDa

Preparation

This antibody was obtained by immunization of HOXD13 recombinant protein (Accession Number: NM_000523). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IF: 1:50-1:500

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:111439(HOXD13 Antibody) at dilution of 1:300



Immunofluorescent analysis of (10% Formaldehyde) fixed A549 cells using Catalog No:111439(HOXD13 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)