

HOXD4 antibody



Catalog Number: 111441

Product name

HOXD4 antibody

Immunogen

Specificity

Human; other species not tested.

Antibody description

HOXD4 Rabbit Polyclonal antibody. Positive WB detected in A549 cells. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of HOXD4 recombinant protein (Accession Number: NM_014621). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:111441(HOXD4 antibody) at dilution of 1:1000