

HP1BP3 antibody



Catalog Number: 111444

Product name

HP1BP3 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse; other species not tested.

Applications

ELISA, WB

Antibody description

HP1BP3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, MCF-7 cells, mouse thymus tissue. Observed molecular weight by Western-blot: 70-75 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of HP1BP3 recombinant protein (Accession Number: BC045660). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111444(HP1BP3 Antibody) at dilution of 1:600

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality