

HPCAL4 antibody

Catalog Number: 111445

Product name

HPCAL4 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

HPCAL4 Rabbit Polyclonal antibody. Positive IHC detected in human medulloblastoma. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 22kd

Preparation

This antibody was obtained by immunization of HPCAL4 recombinant protein (Accession Number: NM_016257). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111445(HPCAL4 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human medulloblastoma using Catalog No:111445(HPCAL4 antibody) at dilution of 1:100 (under 10x lens)