

# Hsc70 antibody



Catalog Number: 111472

## Product name

Hsc70 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

Hsc70 Rabbit Polyclonal antibody. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells, A431 cells, C6 cells, COLO 320 cells, HeLa cells, HepG2 cells, MCF7 cells, NIH/3T3 cells, rat brain tissue. Positive IHC detected in human ovary tumor tissue, human gliomas tissue. Observed molecular weight by Western-blot: 70-73 kDa

## Preparation

This antibody was obtained by immunization of Hsc70 recombinant protein (Accession Number: BC007276). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, IP, WB

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


IP: 1:200-1:2000

IHC: 1:20-1:200

## Validations



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:111472(HSPA8 antibody) at dilution of 1:100



(under 25x lens)

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111472(Hsc70 Antibody) at dilution of 1:2000

# Hsc70 antibody

Catalog Number: 111472

---

IP Result of anti-Hsc70 (IP:Catalog No:111472, 4ug; Detection:Catalog No:111472 1:800) with HEK-293 cells lysate 3680ug.

