

# HMP19 antibody



Catalog Number: 111493

## Product name

HMP19 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

HMP19 Rabbit Polyclonal antibody. Positive IHC detected in human hepat°Cirrhosis. Positive IF detected in HepG2.

## Preparation

This antibody was obtained by immunization of HMP19 recombinant protein (Accession Number: NM\_015980). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, IF


## Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100


## Validations



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:111493(HMP19 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:111493(HMP19 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using HMP19 antibody Catalog No:111493 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).