

HMP19 antibody



Catalog Number: 111493

Product name

HMP19 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

HMP19 Rabbit Polyclonal antibody. Positive IHC detected in human hepat°Cirrhosis. Positive IF detected in HepG2.

Preparation

This antibody was obtained by immunization of HMP19 recombinant protein (Accession Number: NM_015980). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:111493(HMP19 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human hepatocirrhosis using Catalog No:111493(HMP19 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using HMP19 antibody Catalog No:111493 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).