

# HN1 antibody



Catalog Number: 111494

## Product name

HN1 antibody

WB: 1:200-1:2000

IP: 1:500-1:5000

## Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

## Antibody description

HN1 Rabbit Polyclonal antibody. Positive IP detected in MCF-7 cells. Positive WB detected in HeLa cells, A431 cells, HepG2 cells, human brain tissue, human testis tissue, LNCaP cells, mouse colon tissue, mouse testis tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 23-30 kDa

## Validations



## Preparation

This antibody was obtained by immunization of HN1 recombinant protein (Accession Number: NM\_016185). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:111494(HN1 antibody) at dilution of 1:1000

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


ELISA, WB, IF, IP

## Dilutions

Recommended Dilution:



IP Result of anti-HN1 (IP:Catalog No:111494, 4ug; Detection:Catalog No:111494 1:1000) with MCF-7 cells lysate 3200ug.



# HN1 antibody



Catalog Number: 111494

---

Immunofluorescent analysis of HepG2 cells using  
Catalog No:111494(HN1 Antibody) at dilution of

1:25 and Alexa Fluor 488-conjugated AffiniPure  
Goat Anti-Rabbit IgG(H+L)