

HNF1B antibody



Catalog Number: 111496

Product name

HNF1B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HNF1B Rabbit Polyclonal antibody. Positive WB detected in L02 cells, Jurkat cells. Positive IHC detected in human pancreas cancer tissue, human liver cancer tissue. Observed molecular weight by Western-blot: 58-60kd

Preparation

This antibody was obtained by immunization of HNF1B recombinant protein (Accession Number: NM_000458). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:111496(HNF1B antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111496(HNF1B antibody) at dilution of 1:50 (under 10x lens)