

HNF1B antibody

Catalog Number: 111496

Product name

HNF1B antibody

ELISA, WB, IHC

Specificity

Human, Mouse, Rat; other species not tested.

Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Antibody description

HNF1B Rabbit Polyclonal antibody. Positive WB detected in L02 cells, Jurkat cells. Positive IHC detected in human pancreas cancer tissue, human liver cancer tissue. Observed molecular weight by Western-blot: 58-60kd

Validations



Preparation

This antibody was obtained by immunization of HNF1B recombinant protein (Accession Number: NM_000458). Purification method: Antigen affinity purified.

L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:111496(HNF1B antibody) at dilution of 1:800

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111496(HNF1B antibody) at dilution of 1:50 (under 10x lens)