

# HNRNPA2B1 antibody



Catalog Number: 111501

## Product name

HNRNPA2B1 antibody

Recommended Dilution:

WB: 1:500-1:5000

## Specificity

Human, Mouse, Rat; other species not tested.

IP: 1:200-1:2000


IHC: 1:20-1:200

## Antibody description

HNRNPA2B1 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, A549 cells, HeLa cells, Jurkat cells, K-562 cells, Raji cells. Positive IP detected in mouse liver tissue. Positive IHC detected in human lymphoma tissue, human cervical cancer tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 36 kDa and 38 kDa

IF: 1:50-1:500

## Validations



## Preparation

This antibody was obtained by immunization of HNRNPA2B1 recombinant protein (Accession Number: BC000506). Purification method: Antigen affinity purified.

HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:111501(HNRNPA2B1 antibody) at dilution of

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal


## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP, IF

## Dilutions




1:400

IP Result of anti-HNRNPA2B1 (IP:Catalog No:111501, 4ug; Detection:Catalog No:111501 1:500) with mouse liver tissue lysate 4000ug.


# HNRNPA2B1 antibody




Catalog Number: 111501



Immunohistochemistry of paraffin-embedded human lymphoma tissue slide using Catalog No:111501(HNRNPA2B1 Antibody) at dilution of 1:50 (under 40x lens)



Immunohistochemistry of paraffin-embedded human cervical cancer tissue slide using Catalog No:111501(HNRNPA2B1 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog No:111501(HNRNPA2B1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)