

HNRNPUL1 antibody



Catalog Number: 111514

Product name

HNRNPUL1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HNRNPUL1 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in A549 cells, HEK-293 cells, HeLa cells, MCF7 cells, mouse brain tissue. Positive FC detected in HepG2 cells. Positive IF detected in HepG2 cells. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 90-115 kDa

Preparation

This antibody was obtained by immunization of HNRNPUL1 recombinant protein (Accession Number: NM_007040). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF, FC, IP

Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:1000-1:10000


IHC: 1:20-1:200

IF: 1:20-1:200

Validations




A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:111514(HNRNPUL1 antibody) at dilution of



1:2000
Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:111514(HNRNPUL1 antibody) at dilution of 1:100 (under 25x lens)


HNRNPUL1 antibody

Catalog Number: 111514




Immunofluorescent analysis of HepG2 cells, using HNRNPUL1 antibody Catalog No:111514 at 1:50 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent

DNA dye).



IP Result of anti-HNRNPUL1 (IP:Catalog No:111514, 5ug; Detection:Catalog No:111514 1:2000) with HeLa cells lysate 2000ug.



1X10⁶ HepG2 cells were stained with 0.2ug HNRNPUL1 antibody (Catalog No:111514, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.