

HOMER3 antibody



Catalog Number: 111519

Product name

HOMER3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

HOMER3 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 40kd, 45kd

Preparation

This antibody was obtained by immunization of HOMER3 recombinant protein (Accession Number: XM_047439733). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111519(HOMER3 antibody) at dilution of 1:500