

HOMEZ antibody



Catalog Number: 111520

Product name

HOMEZ antibody

Specificity

Human, Rat; other species not tested.

Antibody description

HOMEZ Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, rat kidney tissue. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of HOMEZ recombinant protein (Accession Number: BC130392). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111520(HOMEZ Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human testis slide using Catalog No:111520(HOMEZ Antibody) at dilution of 1:50