

# HOXA2 antibody



Catalog Number: 111528

## Product name

HOXA2 antibody

Polyclonal

## Ig Type

## Specificity

Human, Rat; other species not tested.

Rabbit IgG

## Applications

## Antibody description

HOXA2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, C6 cells, Jurkat cells. Observed molecular weight by Western-blot: 41-43 kDa

ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Preparation

This antibody was obtained by immunization of HOXA2 recombinant protein (Accession Number: NM\_006735). Purification method: Antigen affinity purified.

## Formulation


PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:111528(HOXA2 Antibody) at dilution of 1:600