

HOXA2 antibody



Catalog Number: 111528

Product name

HOXA2 antibody

Specificity

Human, Rat; other species not tested.

Antibody description

HOXA2 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, C6 cells, Jurkat cells. Observed molecular weight by Western-blot: 41-43 kDa

Preparation

This antibody was obtained by immunization of HOXA2 recombinant protein (Accession Number: NM_006735). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:111528(HOXA2 Antibody) at dilution of 1:600