

# HOXA9 antibody

Catalog Number: 111530

## Product name

HOXA9 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

HOXA9 Rabbit Polyclonal antibody. Positive IP detected in mouse appendix tissue. Positive WB detected in mouse colon tissue. Observed molecular weight by Western-blot: 40 kDa

## Preparation

This antibody was obtained by immunization of HOXA9 recombinant protein (Accession Number: NM\_152739). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:111530(HOXA9 Antibody) at dilution of 1:600



IP Result of anti-HOXA9 (IP:Catalog No:111530, 3ug; Detection:Catalog No:111530 1:300) with mouse appendix tissue lysate 4000ug.