

HOXC10 antibody



Catalog Number: 111535

Product name

HOXC10 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

HOXC10 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, A549 cells. Observed molecular weight by Western-blot: 50kd

Preparation

This antibody was obtained by immunization of HOXC10 recombinant protein (Accession Number: NM_017409). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111535(HOXC10 antibody) at dilution of 1:500