

HOXC8 antibody



Catalog Number: 111539

Product name

HOXC8 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HOXC8 Rabbit Polyclonal antibody. Positive WB detected in COLO 320 cells, MDA-MB-453s cells. Observed molecular weight by Western-blot: 33 kDa

Preparation

This antibody was obtained by immunization of HOXC8 recombinant protein (Accession Number: NM_022658). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:111539(HOXC8 Antibody) at dilution of 1:600