

HSD11B1 antibody



Catalog Number: 111542

Product name

HSD11B1 antibody

Polyclonal

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat; other species not tested.

Applications

ELISA, WB

Antibody description

HSD11B1 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, Jurkat cells, mouse liver tissue. Observed molecular weight by Western-blot: 65-70 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Preparation

This antibody was obtained by immunization of HSD11B1 recombinant protein (Accession Number: NM_005525). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:111542(HSD11B1 antibody) at dilution of 1:400

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality