

HSF4 antibody



Catalog Number: 111559

Product name

HSF4 antibody

Specificity

human, mouse; other species not tested.

Antibody description

HSF4 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, human skeletal muscle tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 53-58kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: BC130383). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111559(HSF4 antibody) at dilution of 1:300