

HSFY1 antibody



Catalog Number: 111562

Product name

HSFY1 antibody

Specificity

Human; other species not tested.

Antibody description

HSFY1 Rabbit Polyclonal antibody. Positive WB detected in human testis tissue. Positive IHC detected in human testis tissue. Observed molecular weight by Western-blot: 45 kDa

Preparation

This antibody was obtained by immunization of HSFY1 recombinant protein (Accession Number: XM_017030031). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:111562(HSFY1 Antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:111562(HSFY1 Antibody) at dilution of 1:50 (under 40x lens)