

HSP20 antibody



Catalog Number: 111564

Product name

HSP20 antibody

Specificity

Human, Rat; other species not tested.

Antibody description

HSP20 Rabbit Polyclonal antibody. Positive WB detected in rat heart tissue. Observed molecular weight by Western-blot: 18-20 kDa

Preparation

This antibody was obtained by immunization of HSP20 recombinant protein (Accession Number: NM_144617). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



rat heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111564(HSPB6 Antibody) at dilution of 1:1000