

HYAL3 antibody



Catalog Number: 111577

Product name

HYAL3 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

Immunogen

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HYAL3 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human prostate cancer tissue. Positive WB detected in human brain tissue, HepG2 cells. Observed molecular weight by Western-blot: 47 kDa, 55 kDa

Preparation

This antibody was obtained by immunization of HYAL3 recombinant protein (Accession Number: NM_001200029). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111577(HYAL3 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human testis using Catalog No:111577(HYAL3 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human testis using Catalog No:111577(HYAL3

HYAL3 antibody



Catalog Number: 111577

antibody) at dilution of 1:50 (under 40x lens)