

ICA1L antibody



Catalog Number: 111584

Product name

ICA1L antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

human, mouse, rat; other species not tested.

Antibody description

ICA1L Rabbit Polyclonal antibody. Positive IF detected in Hela. Positive IHC detected in human brain. Positive WB detected in mouse brain tissue, mouse pancreas tissue. Observed molecular weight by Western-blot: 69kd

Preparation

This antibody was obtained by immunization of ICA1L recombinant protein (Accession Number: NM_138468). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111584(ICA1L antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human brain using Catalog No:111584(ICA1L antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human brain using Catalog No:111584(ICA1L

ICA1L antibody

Catalog Number: 111584

antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells, using ICA1L antibody Catalog No:111584 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).