

# ICOS antibody



Catalog Number: 111595

## Product name

Polyclonal

ICOS antibody

## Ig Type

## Specificity

Rabbit IgG

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

ICOS Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, HL-60 cells. Observed molecular weight by Western-blot: 28 kDa, 45-50 kDa

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Preparation

This antibody was obtained by immunization of ICOS recombinant protein (Accession Number: NM\_012092). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111595(ICOS Antibody) at dilution of 1:300

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality