

# ICOSLG antibody



Catalog Number: 111596

## Product name

Polyclonal

ICOSLG antibody

## Ig Type

## Specificity

Rabbit IgG

human; other species not tested.

## Applications

## Antibody description

ELISA, WB

ICOSLG Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 50kd,60kd,70kd

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of ICOSLG recombinant protein (Accession Number: NM\_015259). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111596(ICOSLG antibody) at dilution of 1:1500