

IDI1 antibody

Catalog Number: 111607

Product name

IDI1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

IDI1 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, human liver tissue, L02 cells. Positive IP detected in mouse liver tissue. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of IDI1 recombinant protein (Accession Number: BC019227). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:111607(IDI1 antibody) at dilution of 1:300



IP Result of anti-IDI1 (IP:Catalog No:111607, 4ug; Detection:Catalog No:111607 1:300) with mouse liver tissue lysate 4000ug.