

IER2 antibody



Catalog Number: 111613

Product name

IER2 antibody

Specificity

Human; other species not tested.

Antibody description

IER2 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of IER2 recombinant protein (Accession Number: NM_004907). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111613(IER2 Antibody) at dilution of 1:700