

# IFNAR1 antibody



Catalog Number: 111626

## Product name

IFNAR1 antibody

## Specificity

Human; other species not tested.

## Antibody description

IFNAR1 Rabbit Polyclonal antibody. Positive FC detected in K-562 cells. Positive WB detected in A431 cells, human liver tissue. Observed molecular weight by Western-blot: 64kd

## Preparation

This antibody was obtained by immunization of IFNAR1 recombinant protein (Accession Number: NM\_000629). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, FC

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:111626(IFNAR1 antibody) at dilution of 1:600



1X10<sup>6</sup> K-562 cells were stained with 0.2ug IFNAR1 antibody (Catalog No:111626, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.