

# IGSF2 antibody



Catalog Number: 111634

## Product name

IGSF2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

IGSF2 Rabbit Polyclonal antibody. Positive FC detected in RAW 264.7 cells. Positive WB detected in RAW 264.7 cells. Observed molecular weight by Western-blot: 115 kDa

## Preparation

This antibody was obtained by immunization of IGSF2 recombinant protein (Accession Number: NM\_001256106). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, FC

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Validations



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:111634(IGSF2 Antibody) at dilution of 1:300



1X10<sup>6</sup> RAW 264.7 cells were stained with 0.2ug IGSF2 antibody (Catalog No:111634, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.