

IKZF5 antibody



Catalog Number: 111651

Product name

IKZF5 antibody

Specificity

Human; other species not tested.

Antibody description

IKZF5 Rabbit Polyclonal antibody. Positive IP detected in Jurkat cells. Positive WB detected in Jurkat cells, HeLa cells, HL-60 cells. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of IKZF5 recombinant protein (Accession Number: BC022564). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:111651(IKZF5 Antibody) at dilution of 1:300



IP Result of anti-IKZF5 (IP:Catalog No:111651, 4ug; Detection:Catalog No:111651 1:300) with Jurkat cells lysate 2800ug.