

IL11 antibody



Catalog Number: 111658

Product name

IL11 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

IL11 Rabbit Polyclonal antibody. Positive WB detected in human plasma tissue. Observed molecular weight by Western-blot: 25 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000641). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



human plasma tissue were subjected to SDS PAGE followed by western blot with Catalog No:111658(IL11 Antibody) at dilution of 1:300