

# IL11RA antibody

Catalog Number: 111659

## Product name

IL11RA antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

IL11RA Rabbit Polyclonal antibody. Positive IF detected in N2a cells. Positive FC detected in K-562 cells.

## Preparation

This antibody was obtained by immunization of IL11RA recombinant protein (Accession Number: NM\_001142784). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IF, FC


## Dilutions

Recommended Dilution:


IF: 1:20-1:200

## Validations

N2a neuroblastoma cells fixed under native (4% PFA- Triton X 100) or denaturing (4% PFA + Guanidine thiocyanate) conditions. Labeled with Anti-IL11ra (Protein tech #10264-1-AP)



IF result of anti-IL11RA (Catalog No:111659, 1:200) with N2a cells (4%PFA fixed) by Miss Masue Marbiah, PhD student MRC Prion Unit, Department of Neurodegenerative Diseases, UCL Institute of Neurology.



1X10<sup>6</sup> K-562 cells were stained with 0.2ug IL11RA antibody (Catalog No:111659, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.