

IFT57 antibody

Catalog Number: 111670

Product name

IFT57 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

IFT57 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue, human pancreas tissue. Positive WB detected in mouse retina tissue, HEK-293 cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 49kd

Preparation

This antibody was obtained by immunization of IFT57 recombinant protein (Accession Number: NM_018010). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse retina tissue were subjected to SDS PAGE followed by western blot with Catalog No:111670(IFT57 antibody) at dilution of 1:1000



IP Result of anti-IFT57 (IP:Catalog No:111670, 5ug; Detection:Catalog No:111670 1:800) with HEK-293 cells lysate 2000ug.



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog

IFT57 antibody



Catalog Number: 111670

No:111670(IFT57 antibody) at dilution of 1:50 (under 10x lens)