

# IGHMBP2 antibody



Catalog Number: 111697

## Product name

IGHMBP2 antibody

## Specificity

Human; other species not tested.

## Antibody description

IGHMBP2 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells, human brain tissue. Positive IHC detected in human liver cancer tissue. Observed molecular weight by Western-blot: 109 kDa

## Preparation

This antibody was obtained by immunization of IGHMBP2 recombinant protein (Accession Number: NM\_002180). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:111697(IGHMBP2 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog No:111697(IGHMBP2 Antibody) at dilution of 1:50 (under 40x lens)