

HTR6 antibody

Catalog Number: 111713

Product name

HTR6 antibody

Specificity

Human; other species not tested.

Antibody description

HTR6 Rabbit Polyclonal antibody. Positive IHC detected in human brain tissue.

Preparation

This antibody was obtained by immunization of HTR6 recombinant protein (Accession Number: NM_000871). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:111713(HTR6 Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:111713(HTR6 Antibody) at dilution of 1:100 (under 40x lens).