

HURP antibody

Catalog Number: 111726

Product name

HURP antibody

Specificity

Human; other species not tested.

Antibody description

HURP Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells. Observed molecular weight by Western-blot: 95kd

Preparation

This antibody was obtained by immunization of HURP recombinant protein (Accession Number: XM_017021840). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111726(DLGAP5 antibody) at dilution of 1:500