

# HUS1 antibody



Catalog Number: 111727

## Product name

HUS1 antibody

## Specificity

Human; other species not tested.

## Antibody description

HUS1 Rabbit Polyclonal antibody. Positive WB detected in BxPC-3 cells. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of HUS1 recombinant protein (Accession Number: NM\_004507). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC, WB

## Dilutions


Recommended Dilution:

WB: N/AIHC: 1:20-1:200

## Validations



BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:111727(HUS1 antibody) at dilution of 1:0



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:111727(HUS1 antibody) at dilution of 1:50 (under 10x lens)