

IL2RB antibody



Catalog Number: 111733

Product name

IL2RB antibody

Specificity

human; other species not tested.

Antibody description

IL2RB Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 70-75kd

Preparation

This antibody was obtained by immunization of IL2RB recombinant protein (Accession Number: NM_001346222). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111733(IL2RB antibody) at dilution of 1:100