

# IL3 antibody

Catalog Number: 111735

## Product name

IL3 antibody

## Specificity

Human; other species not tested.

## Antibody description

IL3 Rabbit Polyclonal antibody. Positive WB detected in Recombinant protein, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 20 kDa

## Preparation

This antibody was obtained by immunization of IL3 recombinant protein (Accession Number: NM\_000588). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

## Validations



Recombinant protein were subjected to SDS PAGE followed by western blot with Catalog No:111735(IL3 antibody) at dilution of 1:500



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111735(IL3 Antibody) at dilution of 1:1000