

# IL6R antibody

Catalog Number: 111749

## Product name

IL6R antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

IL6R Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue. Positive IP detected in mouse spleen tissue. Observed molecular weight by Western-blot: 80 kDa

## Preparation

This antibody was obtained by immunization of IL6R recombinant protein (Accession Number: NM\_000565). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:111749(IL6R antibody) at dilution of 1:1000



IP Result of anti-IL6R (IP:Catalog No:111749, 4ug; Detection:Catalog No:111749 1:500) with mouse spleen tissue lysate 4000ug.