

ILKAP antibody



Catalog Number: 111754

Product name

ILKAP antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ILKAP Rabbit Polyclonal antibody. Positive WB detected in human heart tissue, HepG2 cells, K-562 cells. Observed molecular weight by Western-blot: 43-45 kDa

Preparation

This antibody was obtained by immunization of ILKAP recombinant protein (Accession Number: NM_030768). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:111754(ILKAP antibody) at dilution of 1:1500