

IL17RC antibody



Catalog Number: 111766

Product name

IL17RC antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

IL17RC Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 45-50kd

Preparation

This antibody was obtained by immunization of IL17RC recombinant protein (Accession Number: BC006411). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111766(IL17RC antibody) at dilution of 1:200