

# IL1RL2 antibody



Catalog Number: 111774

## Product name

IL1RL2 antibody

## Specificity

Human; other species not tested.

## Antibody description

IL1RL2 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells. Positive FC detected in K-562 cells. Observed molecular weight by Western-blot: 65 kDa

## Preparation

This antibody was obtained by immunization of IL1RL2 recombinant protein (Accession Number: NM\_001351447). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, FC

## Dilutions


Recommended Dilution:

WB: 1:200-1:1000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:111774(IL1RL2 Antibody) at dilution of 1:300



1X10<sup>6</sup> K-562 cells were stained with 0.2ug IL1RL2 antibody (Catalog No:111774, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.