

# IL22RA2 antibody



Catalog Number: 111780

## Product name

IL22RA2 antibody

## Specificity

Human; other species not tested.

## Antibody description

IL22RA2 Rabbit Polyclonal antibody. Positive FC detected in HeLa cells. Positive WB detected in MDA-MB-453s cells. Observed molecular weight by Western-blot: 27 kDa

## Preparation

This antibody was obtained by immunization of IL22RA2 recombinant protein (Accession Number: NM\_181309). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, FC

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:111780(IL22RA2 Antibody) at dilution of 1:600



1X10<sup>6</sup> HeLa cells were stained with 0.2ug IL22RA2 antibody (Catalog No:111780, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.