

# IMPAD1 antibody



Catalog Number: 111783

## Product name

IMPAD1 antibody

## Immunogen

[Human IMPAD1 Recombinant protein \(His tag\)](#)

## Specificity

Human; other species not tested.

## Antibody description

IMPAD1 Rabbit Polyclonal antibody. Positive WB detected in 3T3-L1 cells. Observed molecular weight by Western-blot: 39 kDa

## Preparation

This antibody was obtained by immunization of IMPAD1 recombinant protein (Accession Number: NM\_017813). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



3T3-L1 cells were subjected to SDS PAGE followed by western blot with Catalog No:111783(IMPAD1 Antibody) at dilution of 1:300