

# ING4 antibody

Catalog Number: 111795

## Product name

ING4 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

ING4 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells. Positive IP detected in HEK-293 cells. Observed molecular weight by Western-blot: 29kd

## Preparation

This antibody was obtained by immunization of ING4 recombinant protein (Accession Number: NM\_001127582). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:111795(ING4 antibody) at dilution of 1:200



IP Result of anti-ING4 (IP:Catalog No:111795, 3ug; Detection:Catalog No:111795 1:200) with HEK-293 cells lysate 2700ug.