

INPP5B antibody



Catalog Number: 111808

Product name

INPP5B antibody

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

INPP5B Rabbit Polyclonal antibody. Positive WB detected in mouse spleen tissue, HeLa cells, human placenta tissue. Positive IHC detected in human lymphoma tissue. Observed molecular weight by Western-blot: 75 kDa

Preparation

This antibody was obtained by immunization of INPP5B recombinant protein (Accession Number: BC058932). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:111808(INPP5B antibody) at dilution of 1:200



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:111808(INPP5B antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lymphoma using Catalog

INPP5B antibody



Catalog Number: 111808

No:111808(INPP5B antibody) at dilution of 1:50 (under 40x lens)